

# Lamson

## PTS 4000


# Pneumatic Tube System

The high speed transfer system for movement of goods and samples from point A to point B without wasting time



- ◆ Cash
- ◆ Small Parts
- ◆ Documents
- ◆ Blood Specimens
- ◆ Laboratory samples
- ◆ Packing slips, keys, tools, cheques, mail .etc.

**LAMSON**



ESTABLISHED 1898

[www.lamson.com.au](http://www.lamson.com.au)

## Problems and Challenges:

Businesses are always under pressure to streamline business processes, improve efficiencies and add value to the operations through cost cutting and innovative approaches.

## Solution:

Lamson PTS is an innovative solution that not only improves efficiency by cutting costs associated with the manual delivery, but also adds value to your business by working tirelessly for you 24/7.

The Lamson PTS is one worker that does not demand overtime, premium for working unreasonable hours and fast pace. No longer will you be subjected to late deliveries, incorrect deliveries, loss in efficiency, through staff not being available or away on sick leave. Employ Lamson PTS, the ever reliable, low maintenance solutions that will eliminate all inefficiency associated with in-house transportation.



### Why Invest in Lamson PTS

Lamson PTS automates the transport of materials, pathology samples, pharmaceuticals, cash, other valuable items and light goods

This is a low capital cost with multiple user benefits. The return on investment is excellent and will meet the most stringent of evaluations

The very latest opto electronics allow the highest protection from static discharge ensuring trouble free operation.


The unique LCD touch screen is the very latest technology allowing simple utilization for all staff.

Remote diagnostics, e-mail alerts and effortless integration with popular operating system.

### Send Stations

There are several different variations of stations available to suit your needs-

- ◆ Front Loading
- ◆ Top Entry
- ◆ Bottom Entry
- ◆ Through Stations.





During an operation, removed tissue samples can be sent to pathology to be analysed and the results returned to the operating theatre to determine the continuation of the operation with minimal delays to patient and medical staff.

- ◆ Swivel type carriers open at both ends for easy access
- ◆ Carriers are available with electronic transponders that track exactly where the carrier is
- ◆ Carriers can be sent with a request or order for cash, drugs, etc
- ◆ Carriers can be sent into secure cabinets or vacuum safes




- ◆ Cash (coins & notes) can be quickly and securely transported in either direction, as well as receipts and change requests
- ◆ Audit trails are available for staff accountability

#### CARRIER DIAMETERS

- ◆ 75mm / 3 inches
- ◆ 90mm / 3.5 inches
- ◆ 110mm / 4 inches
- ◆ 160mm / 6 inches
- ◆ Other sizes available on request



# MULTI-FUNCTIONAL PNEUMATIC TUBE SYSTEM



Lamson were established in 1898 and ever since then they have moved along with technology and improved on the ways materials and products can be transferred via pneumatic tube system. The latest being PTS.



Lamson Concepts Pty Limited

HEAD OFFICE \_\_\_\_\_

7/14-38 Bellona Avenue,  
Regents Park NSW 2143  
AUSTRALIA

Phone: 61 2 97437322

Fax: 61 2 9743 8569

Email: info@lamson.com.au

OFFICES \_\_\_\_\_

## AUSTRALIA

- Victoria
- Queensland
- South Australia
- Western Australia
- Australian Capital Territory

## NEW ZEALAND

- Auckland

DISTRIBUTORS \_\_\_\_\_

- China
- United States
- Canada
- Indonesia
- Saudi Arabia
- Korea
- South Africa
- Malaysia
- Egypt

[www.lamson.com.au](http://www.lamson.com.au)